



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **GLUCOSE OXIDASE** MANUFACTURE DATE: 12/30/2014
 CATALOG #: G519-1L
 LOT #: 58636 EXP. DATE: 03/2017

<u>COMPONENTS</u>	<u>CATALOG #</u>	<u>LOT #</u>	<u>EXP. DATE</u>
GLUCOSE OXIDASE REAGENT	G519-A	53562	03/2017
GLUCOSE STANDARD (100 mg/dL)	G518-B	55067	05/2017

TOLERANCE LEVELS OF CONTROLS


<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (on TC-Matrix)	<u>TARGET</u> (mg/dL)
27621	12/2015	83.2 – 101.7	92.5
27602	03/2015	245.9 – 300.6	273.3

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27621	92.0	PASS
27602	277.7	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY:


 Jacqueline Sanchez
 Quality Control

DATE:

12/30/2014