

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.  
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880  
Fax: (714) 463-1169  
[www.tecodiag.com](http://www.tecodiag.com)

**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **GLUCOSE OXIDASE** MANUFACTURE DATE: 04/15/2015  
CATALOG #: G519-1L  
LOT #: 60531 EXP. DATE: 01/2018

---

<u>COMPONENTS</u>	<u>CATALOG #</u>	<u>LOT #</u>	<u>EXP. DATE</u>
GLUCOSE OXIDASE REAGENT	G519-A	59036	01/2018
GLUCOSE STANDARD (100 mg/dL)	G518-B	60155	03/2018

---

**TOLERANCE LEVELS OF CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mg/dL)
27621	12/2015	83.2 - 101.7	92.5
27622	12/2015	258.1 - 315.5	286.8

---


**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27621	95.5	PASS
27622	295.2	PASS

---

COMMENTS: VISUAL CHECK OK.

APPROVED BY:

  
Jacqueline Sanchez  
Quality Control

DATE:

04/16/2015