

**TECO DIAGNOSTICS**

1268 N. LAKEVIEW AVE. ANAHEIM, CA 92807, U.S.A.  
PHONE: (714) 463-1122 or (800) 222-9880  
FAX: (714) 463-1169  
[www.tecodiag.com](http://www.tecodiag.com)

**CERTIFICATE OF ANALYSIS**  
**(BATCH RELEASE)**

1. PRODUCT NAME: **GLYCOHEMOGLOBIN**  
ELEVATED CONTROL LOT#: 013002 EXP. DATE: 01/2013
2. METHOD/PROCEDURE: CATION EXCHANGE RESIN
3. DESCRIPTION: LYOPHILIZED CONTROL
4. KEY INGREDIENTS: PREPARED FROM HUMAN ERYTHROCYTES
5. CONTROL MANUFACTURER: TECO
6. MANUFACTURE DATE: 09/14/2010


**TOLERANCE LEVELS OF CONTROLS**

<u>CONTROL OR STANDARD</u>	<u>LOT #</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u>
Standard	052181	10/2012	8.5 – 11.5	10.0
Elevated	013002	01/2013	11.9 – 16.1	14.0

**TEST VALUES/RESULTS**

<u>CONTROL OR STANDARD</u>	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>MEAN</u>	<u># PASS</u>
Standard	10	10	10	10	10	10	10	6/6
Elevated	12.2	12.3	12.9	14.2	14.1	15.5	13.5	6/6

COMMENTS: VISUAL CHECK OK

APPROVE D BY:	 Quality Control	DATE :	09/15/2010
------------------	--	-----------	------------