

**TECO DIAGNOSTICS**

1268 N. LAKEVIEW AVE. ANAHEIM, CA 92807, U.S.A.

PHONE: (714) 463-1122 or (800) 222-9880

FAX: (714) 463-1169

[www.tecodiag.com](http://www.tecodiag.com)

**CERTIFICATE OF ANALYSIS**  
**(BATCH RELEASE)**

1. PRODUCT NAME: **GLYCOHEMOGLOBIN**  
NORMAL CONTROL LOT#:                      EXP. DATE:  
ELEVATED CONTROL LOT#: 102092 EXP. DATE: 10/2012
2. METHOD/PROCEDURE: CATION EXCHANGE RESIN
3. DESCRIPTION: LYOPHILIZED CONTROL
4. KEY INGREDIENTS: PREPARED FROM HUMAN ERYTHROCYTES
5. CONTROL MANUFACTURER: TECO
6. MANUFACTURE DATE: 12/18/2003

**TOLERANCE LEVELS OF CONTROLS**

<u>CONTROL OR STANDARD</u>	<u>LOT #</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u>
Normal	102091	1/2012	5.9 – 8.1	7.0
Standard	052181	10/2012	8.5 – 11.5	10.0
Elevated	102092	10/2012	11.9 – 16.1	14.0

**TEST VALUES/RESULTS**

<u>CONTROL OR STANDARD</u>	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>MEAN</u>	<u># PASS</u>
Normal	7.8	7.8	6.5	6.6	77.0	7.7	7.2	6/6
Standard	9.3	9.5	10.2	10.7	10.1	10.2	10.0	6/6
Elevated	12.7	12.8	13.6	13.8	12.7	13.6	13.2	6/6

COMMENTS: VISUAL CHECK OK