

**TECO DIAGNOSTICS**

1268 N. LAKEVIEW AVE. ANAHEIM, CA 92807, U.S.A.

PHONE: (714) 463-1122 or (800) 222-9880

FAX: (714) 463-1169

[www.tecodiag.com](http://www.tecodiag.com)

**CERTIFICATE OF ANALYSIS**  
**(BATCH RELEASE)**

1. PRODUCT NAME: **GLYCOHEMOGLOBIN**  
NORMAL CONTROL LOT#: 071382 EXP. DATE: 07/2011  
ELEVATED CONTROL LOT#: EXP. DATE:
2. METHOD/PROCEDURE: CATION EXCHANGE RESIN
3. DESCRIPTION: LYOPHILIZED CONTROL
4. KEY INGREDIENTS: PREPARED FROM HUMAN ERYTHROCYTES
5. CONTROL MANUFACTURER: TECO
6. MANUFACTURE DATE: 04/27/09

**TOLERANCE LEVELS OF CONTROLS**

<u>CONTROL OR STANDARD</u>	<u>LOT #</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u>
Normal	071382	07/2011	5.9 – 8.1	7.0
Standard	052181	05/2011	8.5 – 11.5	10.0
Elevated	081581	08/2011	11.9 – 16.1	14.0

**TEST VALUES/RESULTS**

<u>CONTROL OR STANDARD</u>	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>MEAN</u>	<u># PASS</u>
Normal	7.8	8.0	8.1	7.7	7.8	8.1	7.9	6/6
Standard	9.4	9.7	11.0	9.6	9.6	10.8	10.0	6/6
Elevated	12.4	12.3	11.5	13.3	13.2	12.1	12.5	5/6

COMMENTS: VISUAL CHECK OK