



CERTIFICATE OF ANALYSIS
(BATCH RELEASE)

1. PRODUCT NAME: **GLYCOHEMOGLOBIN** CAT. # LOT# EXP. DATE
 RESIN REAGENT G540-A 101001 10/2012
 LYSING REAGENT G540-B 40901 04/2012
 STANDARD (10%) G540-C 101201 10/2013
 NORMAL CONTROL G541-1 62001 06/2013
 ELEVATED CONTROL G542-1 102092 10/2012
2. METHOD/PROCEDURE: CATION EXCHANGE RESIN
3. DESCRIPTION: LIQUID RESIN AND LYSING
 LYOPHILIZED STANDARD AND CONTROLS
4. KEY INGREDIENTS: RESIN: 8 mg/ml CATION EXCHANGE RESIN
 LYSING: 10 mM POTASSIUM CYANIDE
 STANDARD: GLYCOHEMOGLOBIN
5. CONTROL MANUFACTURER: TECO
6. MANUFACTURE DATE: 11/29/2010


TOLERANCE LEVELS OF CONTROLS

<u>CONTROL OR STANDARD</u>	<u>LOT #</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (%)</u>
Normal	62001	06/2013	5.9 – 8.1	7.0
Standard	101201	10/2013	8.5 – 11.5	10.0
Elevated	102092	10/2012	11.9 – 16.1	14.0

TEST VALUES/RESULTS

<u>CONTROL OR STANDARD</u>	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>MEAN</u>	<u># PASS</u>
Normal	6.7	7.7	6.6	7.8	7.6	7.7	7.3	6/6
Standard	10.2	9.9	10.4	9.5	9.4	9.5	10.0	6/6
Elevated	13.3	15.2	13.9	14.9	13.7	14.3	14.2	6/6

COMMENTS: VISUAL CHECK OK.

APPROVE D BY:	 Quality Control	DATE :	<u>11/30/2010</u>
------------------	--	-----------	-------------------