



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **GLYCOHEMOGLOBIN** MANUFACTURE DATE: 08/19/2011
 CATALOG #: G540-40
 LOT #: 34072 EXP. DATE: 06/2013

<u>COMPONENTS</u>	<u>CATALOG #</u>	<u>LOT#</u>	<u>EXP. DATE</u>
1. RESIN REAGENT	G540-A	062711	06/2013
2. LYSING REAGENT	G540-B	062811	06/2013
3. STANDARD (10%)	G540-C	012211	01/2014

CONTROL TOLERANCE LEVELS

<u>CONTROL OR STANDARD</u>	<u>LOT #</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (%)</u>
Standard	012211	01/2014	8.5 – 11.5	10.0
Normal	62001	06/2013	5.9 – 8.1	7.0
Elevated	52701	05/2013	11.9 – 16.1	14.0

TEST VALUES/RESULTS

<u>CONTROL OR STANDARD</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
Standard	11.0	9.0	10.0	2/2
Normal	6.5	6.6	6.6	2/2
Elevated	13.3	13.3	13.3	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	08/19/2011
--------------	--	-------	------------