



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **HDL – CHOLESTEROL (PEG)**      MANUFACTURE DATE: 10/22/2010  
 CATALOG #: H511-20  
 LOT #: 28290      EXP. DATE: 08/2012

COMPONENTS: HDL CHOLESTEROL REAGENT(PEG)    LOT #: 26459    EXP. DATE: 12/2012  
 HDL CHOLESTEROL STANDARD    LOT #: 24650    EXP. DATE: 08/2012

BLANK: 0.011    AT 505 NM

STANDARD: 50 mg/dL

**TOLERANCE LEVELS OF SERUM CONTROL**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> <u>(mg/dL)</u>
<u>921203</u>	<u>07/2011</u>	<u>26(32)38</u>	<u>32</u>
<u>921203</u>	<u>07/2011</u>	<u>99(116)133</u>	<u>116</u>

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
921203	26.4	29.2	27.8	2/2
921203	122.8	125.3	124.05	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:       DATE: 11/10/2010  
Quality Control