



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **HDL – CHOLESTEROL (PEG)**      MANUFACTURE DATE: 12/01/2010  
 CATALOG #: H511-20  
 LOT #: 29077      EXP. DATE: 04/2013

COMPONENTS: HDL CHOLESTEROL REAGENT(PEG)    LOT #: 29080    EXP. DATE: 04/2013  
 HDL CHOLESTEROL STANDARD    LOT #: 29081    EXP. DATE: 04/2013

BLANK: 0.033    AT 505 NM

STANDARD: 50 mg/dL

**TOLERANCE LEVELS OF SERUM CONTROL**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> <u>(mg/dL)</u>
<u>971071-107</u>	<u>03/2011</u>	<u>23.5(29.4)35.3</u>	<u>29.4</u>
<u>971071-107</u>	<u>03/2011</u>	<u>71.4(89.3)107.2</u>	<u>89.3</u>

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
<u>910701-107</u>	27.7	27.6	27.65	2/2
<u>910701-107</u>	101.9	97.1	99.5	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Casey*      DATE: 12/03/2010  
 Quality Control