



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **HDL – CHOLESTEROL (PEG)** MANUFACTURE DATE: 11/10/2011
 CATALOG #: H511-20
 LOT #: 35889 EXP. DATE: 12/2013

COMPONENTS: HDL CHOLESTEROL REAGENT(PEG) LOT #: 33695 EXP. DATE: 12/2013
 HDL CHOLESTEROL STANDARD LOT #: 33696 EXP. DATE: 12/2013

BLANK: 0.011 AT 505 NM

STANDARD: 50 mg/dL

TOLERANCE LEVELS OF SERUM CONTROL

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mg/dL)
111201-112	04/2013	31 - 47	39
111201-112	04/2013	65 - 87	76

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
111201-112	45.5	40.6	43.1	2/2
111201-112	66.5	81.4	74	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 11/11/2011

Nour Qasqas
 Director of Quality Assurance