



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **HDL – CHOLESTEROL (PEG)** MANUFACTURE DATE: 06/19/2014
 CATALOG #: H511-20
 LOT #: 55721 EXP. DATE: 09/2016

COMPONENTS: HDL CHOLESTEROL REAGENT(PEG) LOT #: 53855 EXP. DATE: 09/2016
 HDL CHOLESTEROL STANDARD LOT #: 53856 EXP. DATE: 09/2016

TOLERANCE LEVELS OF SERUM CONTROL

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mg/dL)
14461	03/2016	66.4 – 99.7	83.0
14462	03/2016	19.6 – 29.4	24.5

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
14461	81.2	PASS
14462	25.9	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY:



 Jacqueline Sanchez
 Quality Control

DATE:

