



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **IRON / TIBC** MANUFACTURE DATE: 05/18/2010
 CATALOG #: I592-100
 LOT #: 25369 EXP. DATE: 04/2012

TOLERANCE LEVELS OF SERUM CONTROLS

SERA ID	EXP. DATE	TOLERANCE RANGE		TARGET (ug/dL)	
		IRON	TIBC	IRON	TIBC
15271	03/2012	127 – 207	233 – 349	167	291
15272	03/2012	35 – 59	118 – 176	47	147

TEST VALUES/RESULTS

SERA ID	IRON		TIBC		MEAN		# PASS
	1	2	1	2	IRON	TIBC	
15271	161	146	296.6	260	153.5	265.1	4/4
15272	52	47	175.6	158.7	265.1	157.2	4/4

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 05/18/2010
 Quality Control