



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **IRON / TIBC** MANUFACTURE DATE: 01/04/2011
 CATALOG #: I592-100
 LOT #: 29140 EXP. DATE: 07/2012

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. IRON BUFFER REAGENT	I592-A	29562	01/2013
2. UIBC BUFFER REAGENT	I592-B	24541	07/2012
3. IRON COLOR REAGENT	I592-C	29564	01/2013
4. IRON STANDARD (500 ug/dL)	I592-D	29565	01/2013

TOLERANCE LEVELS OF SERUM CONTROLS

SERA ID	EXP. DATE	TOLERANCE RANGE		TARGET (ug/dL)	
		IRON	TIBC	IRON	TIBC
15281	09/2012	170-220	233 – 349	193	291
15282	09/2012	40.2-60	160-240	50	199.4

TEST VALUES/RESULTS

SERA ID	IRON		TIBC		MEAN		# PASS
	1	2	1	2	IRON	TIBC	
15281	188	203	274	287	195.5	280.5	4/4
15282	43	50	206	205	46.5	205.5	4/4

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Basen* DATE: 01/06/2011
 Quality Control