



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **IRON / TIBC** MANUFACTURE DATE: 02/23/2015
 CATALOG #: 1592-100
 LOT #: 59582 EXP. DATE: 08/2016

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. IRON BUFFER REAGENT	1592-A	55424	08/2016
2. UIBC BUFFER REAGENT	1592-B	54997	08/2016
3. IRON COLOR REAGENT	1592-C	54998	08/2016
4. IRON STANDARD (500 ug/dL)	1592-D	56624	10/2016

TOLERANCE LEVELS OF SERUM CONTROLS


SERA ID	EXP. DATE	TOLERANCE RANGE		TARGET (ug/dL)	
		IRON	TIBC	IRON	TIBC
27621	12/2015	160 - 240	237 - 355	200	296
27602	03/2015	29 - 67	144 - 216	48	180

TEST VALUES/RESULTS

SERA ID	MEAN		PASS/FAIL
	IRON	TIBC	
27621	190.5	275	PASS
27602	40	166	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY: _____


 Roxana Magdaleno
 Quality Control

DATE: _____

02/24/2015