



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **LDH – L LIQUID** MANUFACTURE DATE: 02/11/2010
 CATALOG #: L533-60
 LOT #: 23387 EXP. DATE: 12/2011

COMPONENTS: 1. LDH-L LIQUID REAGENT 1 LOT#: 19803 EXP. DATE: 12/2011

REAGENT BLANK VS. WATER: 0.357 AT 340 NM

TOLERANCE LEVELS OF SERUM CONTROL

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u>
15271	03/2012	115-173	144
15272	03/2012	273-409	341

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15271	132.6	146.3	139.5	2/2
15272	314.5	347.0	330.8	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Casper* DATE: 02/11/2010
 Quality Control