



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **LDH-L (KINETIC)**                      MANUFACTURE DATE: 01/20/2011  
 CATALOG #: L535-60  
 LOT #: 28398    EXP. DATE: 12/2012

<u>COMPONENT</u>	<u>CATALOG#</u>	<u>LOT#</u>	<u>EXP. DATE</u>
LDH REAGENT	L535-A	26761	12/2012

REAGENT BLANK VS. WATER: 0.331 AT 340 NM

**TOLERANCE LEVELS OF CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (U/L)</u>
15271	03/20012	115-173	144
15272	03/2012	269-403	336

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15271	116.0	119.9	117.95	2/2
15272	304.2	335.3	319.75	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED  
BY:

  
Quality Control

DATE: 01/20/2011