



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **LDH-L (KINETIC)**                      MANUFACTURE DATE: 01/11/2011  
 CATALOG #: L535-60  
 LOT #: 30034    EXP. DATE: 04/2013

<u>COMPONENT</u>	<u>CATALOG#</u>	<u>LOT#</u>	<u>EXP. DATE</u>
LDH REAGENT	L535-A	29169	04/2013

REAGENT BLANK VS. WATER: 0.335 AT 340 NM

**TOLERANCE LEVELS OF CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (U/L)</u>
15271	03/20012	113.3-170	141.7
15272	03/2012	279-418.4	348.7

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15271	163	162	162.5	2/2
15272	384	388	386	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:                       DATE: 01/11/2011