



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **LDH-L (KINETIC)** MANUFACTURE DATE: 08/28/2014
 CATALOG #: L535-60
 LOT #: 56855 EXP. DATE: 01/2017

<u>COMPONENT</u>	<u>CATALOG#</u>	<u>LOT#</u>	<u>EXP. DATE</u>
LDH REAGENT	L535-A	56314	01/2017

TOLERANCE LEVELS OF CONTROLS


<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-Matrix)	<u>TARGET</u> (U/L)
27611	05/2015	129.5 – 194.3	161.9
27602	03/2015	307.3 – 460.9	384.1

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27611	175.0	PASS
27602	407.7	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY:


 Jacqueline Sanchez
 Quality Control

DATE:

8/28/2014