



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **LDH-L (KINETIC)** MANUFACTURE DATE: 01/21/2015
 CATALOG #: L535-60
 LOT #: 58980 EXP. DATE: 07/2017

<u>COMPONENT</u>	<u>CATALOG#</u>	<u>LOT#</u>	<u>EXP. DATE</u>
LDH REAGENT	L535-A	58719	07/2017

TOLERANCE LEVELS OF CONTROLS


<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-Matrix)	<u>TARGET</u> (U/L)
27621	12/2015	134.2 – 201.3	167.7
27602	03/2015	307.3 – 460.9	384.1

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27621	168.4	PASS
27602	404.3	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY:


 Jacqueline Sanchez
 Quality Control

DATE:

01/21/2015