



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **INORGANIC PHOSPHORUS REAGENT SET (UV METHOD)** MANUFACTURE DATE: 10/22/2008
 CATALOG #: 1516-480
 LOT #: 17005 EXP. DATE: 01/2011

<u>COMPONENTS:</u>	<u>CAT. #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. INORGANIC PHOSPHORUS UV REAGENT	I516-A	16303	01/2011
2. INORGANIC PHOSPHORUS UV STANDARD (5.0 mg/dL)	I516-B	16304	01/2011

BLANK: 0.084 AT 340 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (mg/dL)</u>
15221	08/2009	2.7 – 3.5	3.1
15222	08/2009	5.0 – 6.2	5.6

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15221	3.0	3.0	3.0	2/2
15222	5.6	5.6	5.6	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *gdkta* DATE: 10/22/2008
 QA/QC