



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **INORGANIC PHOSPHORUS REAGENT SET (UV METHOD)** MANUFACTURE DATE: 5/28/09
 CATALOG #: 1516-480
 LOT #: 19564 EXP. DATE: 11/2011

<u>COMPONENTS:</u>	<u>CAT. #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. INORGANIC PHOSPHORUS UV REAGENT	I516-A	19562	11/2011
2. INORGANIC PHOSPHORUS UV STANDARD (5.0 mg/dL)	I516-B	19563	11/2011

BLANK: .093 AT 340 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mg/dL)
15261	08/2011	2.6-3.4	3.0
15262	08/2011	4.9-6.3	5.6

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15261	3.3	3.3	3.3	2/2
15262	6.1	6.0	6.1	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:


 Quality Control

DATE: 5/28/09