



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **INORGANIC PHOSPHORUS REAGENT SET (UV METHOD)** MANUFACTURE DATE: 02/07/2011
 CATALOG #: 1516-480
 LOT #: 29412 EXP. DATE: 03/2012

<u>COMPONENTS:</u>	<u>CAT. #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. INORGANIC PHOSPHORUS UV REAGENT	I516-A	20627	03/2012
2. INORGANIC PHOSPHORUS UV STANDARD (5.0 mg/dL)	I516-B	20628	03/2012

BLANK: 0.163 AT 340 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mg/dL)
15281	09/2012	2.0-3.0	2.5
15272	03/2012	4.9-6.3	5.6

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	2.4	2.4	2.4	2/2
15272	5.3	5.0	5.15	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 02/07/2011
 BY: Quality Control