



CERTIFICATE OF ANALYSIS

1. PRODUCT NAME: **POTASSIUM**
 REAGENT CATALOG # LOT # EXP. DATE
 P605-50 21061 11/2012

2. METHOD/PROCEDURE: END POINT, TURBIDITY

3. DESCRIPTION: CLEAR LIQUID

4. KEY INGREDIENTS: REAGENT: SODIUM TETRAPHENYL BORON

5. REAGENT BLANK: 0.031 AT 505 NM

6. CONTROL MANUFACTURER: TECO

7. MANUFACTURE DATE: 04/28/2010

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID/LOT #</u>	<u>EXP. DATE</u>	<u>RANGE</u> (On TC-84)	<u>TARGET</u> (mEq/L)
15271	03/2012	3.7 - 4.7	4.3
15272	03/2012	4.9 - 6.9	5.9

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>MEAN</u>	<u># PASS</u>
15271	4.3	4.1	4.4	4.3	4.2	4.4	4.3	6/6
15272	5.8	5.7	5.7	5.8	5.6	5.6	5.7	6/6

COMMENTS: VISUAL CHECK OK.

APPROVED BY:


Quality Control

DATE: 04/28/2010