



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **POTASSIUM** MANUFACTURE DATE: 11/12/2010
CATALOG #: P605-50
LOT #: 29010 EXP. DATE: 04/2013

COMPONENTS: 1. POTASSIUM REAGENT LOT #: D100501022 EXP. DATE: 05/2013
2. POTASSIUM STANDARD LOT #: 29013 EXP. DATE: 04/2013

BLANK: 0.035 AT 505 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mEq/L)
15271	03/2012	3.7 – 4.7	4.3
15272	03/2012	4.9 – 6.9	5.9

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	3.9	3.9	3.9	2/2
15272	5.4	5.2	5.3	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:


Quality Control

DATE: 11/15/2010