



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **POTASSIUM**  
 CATALOG #: P605-50  
 LOT #: 56173

MANUFACTURE DATE: 07/16/2014

EXP. DATE: 09/2016

---

COMPONENTS: 1. POTASSIUM REAGENT LOT #: 54227 EXP. DATE: 10/2016  
 2. POTASSIUM STANDARD LOT #: 52781 EXP. DATE: 09/2016

---

**TOLERANCE LEVELS OF SERUM CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (on TC-Matrix)	<u>TARGET</u> (mmol/l)
27611	05/2015	3.5 – 4.3	3.9
27602	03/2015	5.1 – 7.7	6.4

---

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27611	4.0	PASS
27602	6.3	PASS

---

COMMENTS: VISUAL CHECK OK.

APPROVED BY: \_\_\_\_\_

*Susan Tang*  
 Susan Tang  
 Quality Control

DATE: \_\_\_\_\_

*07/17/2014*