



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **SODIUM** MANUFACTURE DATE: 04/18/2011
 CATALOG #: S600-50
 KIT LOT #: 31991 EXP. DATE: 11/2013

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. SODIUM FILTRATE REAGENT	S600-A	31312	11/2013
2. SODIUM ACID REAGENT	S600-B	31979	11/2013
3. SODIUM COLOR REAGENT	S600-C	31980	11/2013
4. SODIUM STANDARD (150 mEq/L)	S600-D	31981	11/2013

BLANK: 1.522 AT 550 NM


TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mEq/L)
15281	09/2012	88-133	111
15272	03/2012	114 – 144	129

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	101.7	98.6	100.2	2/2
15272	117.3	117.1	117.2	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	04/18/2011
-----------------	--	-------	------------