



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **SODIUM** MANUFACTURE DATE: 07/19/2011
 CATALOG #: S600-50
 KIT LOT #: 33740 EXP. DATE: 12/2013

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. SODIUM FILTRATE REAGENT	S600-A	33523	12/2013
2. SODIUM ACID REAGENT	S600-B	33524	12/2013
3. SODIUM COLOR REAGENT	S600-C	33525	12/2013
4. SODIUM STANDARD (150 mEq/L)	S600-D	33526	12/2013

BLANK: 1.861 AT 550 NM

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mEq/L)
15281	09/2012	88-133	111
15282	09/2012	101-152	126

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	90.7	92.5	91.6	2/2
15282	118.5	118.7	118.6	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	<u>07/19/2011</u>
-----------------	--	-------	-------------------