



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **SODIUM** MANUFACTURE DATE: 11/18/2011
 CATALOG #: S600-50
 KIT LOT #: 36100 EXP. DATE: 04/2014

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. SODIUM FILTRATE REAGENT	S600-A	35737	04/2014
2. SODIUM ACID REAGENT	S600-B	35738	04/2014
3. SODIUM COLOR REAGENT	S600-C	35739	04/2014
4. SODIUM STANDARD (150 mEq/L)	S600-D	35740	04/2014

BLANK: 1.705 AT 550 NM

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mEq/L)
15281	09/2012	88-133	111
15282	09/2012	101 - 152	126

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	109.6	103.1	106.4	2/2
15282	130.9	125.9	128.4	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: _____

Nour Qasqas
 Nour Qasqas
 Director of Quality Assurance

DATE: 11/18/2011