



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **BILIRUBIN TOTAL** MANUFACTURE DATE: 11/18/2009
 CATALOG #: B576-480
 LOT #: 22370 EXP. DATE: 6/2012

<u>COMPONENTS:</u>	<u>CATALOG #</u>	<u>LOT#</u>	<u>EXP DATE</u>
1. TOTAL BILIRUBIN REAGENT	B576-A	21967	11/2012
2. NITRITE REAGENT	B538-B	19781	06/2012
3. CALIBRATOR (5.0 mg/dl)	B538-C	19780	06/2012

BLANK: 0.067 AT 550 nm

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dl)
15271	03/2012	0.4 – 1.9	1.2
15272	03/2012	4.9 – 7.3	6.1

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	0.6	0.9	0.8	2/2
15272	6.3	6.0	6.2	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Pasquini* DATE: 11/18/2009
Quality Control