



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **TOTAL PROTEIN** MANUFACTURE DATE: 04/25/2014
 CATALOG #: T528-480
 LOT #: 54450 EXP. DATE: 08/2015

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT #:</u>	<u>EXP. DATE:</u>
1. TOTAL PROTEIN REAGENT	T528-A	52768	03/2017
2. TOTAL PROTEIN STANDARD (5.0 g/dL)	T528-B	41877	08/2015

TOLERANCE LEVELS OF SERUM CONTROLS


<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (on TC Matrix)	<u>TARGET</u> (g/dL)
27611	05/2015	6.4 – 7.9	7.1
27602	03/2015	4.1 – 5.1	4.6

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27611	7.0	PASS
27602	4.3	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY:



 Jacqueline Sanchez
 Quality Control

DATE:

04/28/2014