



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **TOTAL PROTEIN**                      MANUFACTURE DATE:    03/16/2015  
 CATALOG #:    T528-480  
 LOT #:                      59955                                      EXP. DATE:                      03/2017

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT #:</u>	<u>EXP. DATE:</u>
1. TOTAL PROTEIN REAGENT	T528-A	52768	03/2017
2. TOTAL PROTEIN STANDARD (5.0 g/dL)	T528-B	57373	09/2017

**TOLERANCE LEVELS OF SERUM CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (on TC Matrix)	<u>TARGET</u> (g/dL)
27621	12/2015	6.0 - 7.4	6.7
27602	03/2015	4.1 - 5.1	4.6

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27621	7.2	PASS
27602	4.5	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *Jenifer Ohta*  
 Jenifer Ohta  
 Quality Assurance Manager

DATE: 03/16/2015