



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **TOTAL PROTEIN**                      MANUFACTURE DATE: 10/12/2009  
 CATALOG #: T528-100  
 LOT #: 19559    EXP. DATE: 6/2012

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT #:</u>	<u>EXP. DATE:</u>
1. TOTAL PROTEIN REAGENT	T528-A	19557	6/2012
2. TOTAL PROTEIN STANDARD (5.0 g/dL)	T528-B	21078	10/2012

BLANK: 0.097 AT 550 NM

**TOLERANCE LEVELS OF SERUM CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (g/dL)
15271	03/2012	4.8 – 7.0	5.9
15272	03/2012	3.2 – 4.8	4.0

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	4.9	4.8	4.9	2/2
15272	3.6	3.6	3.6	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Bassey*                      DATE: 10/12/2009  
Quality Control