



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: TC SERUM CONTROL LEVEL I MANUFACTURE DATE: 12/24/2008

LOT #: 17859 EXP. DATE: 08/2011

COMPONENT: TC SERUM CONTROL LEVEL I LOT#: 15261 EXP. DATE: 08/2011

REAGENT (ALBUMIN) BLANK VS WATER: 0.387 AT 600 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u>
15261	08/2011	<u>URIC ACID</u> 3.2 – 4.2	3.7

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15261	3.7	3.8	3.75	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *J. Ohta* DATE: 12/24/2008
 QA/QC