

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **TRIGLYCERIDE-GPO (LIQUID)** MANUFACTURE DATE: 03/18/2015
CATALOG #: T532-480
LOT #: 60071 EXP. DATE: 06/2016

COMPONENTS: 1. REAGENT LOT#: 58332 EXP. DATE: 06/2016
 2. STANDARD LOT#: 56680 EXP. DATE: 08/2017

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> <u>(TC MATRIX)</u>	<u>TARGET</u> <u>(mg/dL)</u>
27621	12/2015	139.4 – 209.1	174.2
27602	03/2015	50.6 – 84.3	67.4

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u># PASS</u>
27621	157.4	2/2
27602	61.0	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: jenifer Ohta
Jenifer Ohta
Quality Assurance Manager

DATE: 03/18/2015