

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: TRIGLYCERIDE-GPO (LIQUID) MANUFACTURE DATE: 07/23/2015
CATALOG #: T532-480
LOT #: 62035 EXP. DATE: 01/2017

COMPONENTS: 1. REAGENT LOT#: 61789 EXP. DATE: 01/2017
2. STANDARD LOT#: 56680 EXP. DATE: 08/2017

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (Manual Method)	<u>TARGET</u> (mg/dL)
27621	12/2015	139.4 – 209.1	174.2
27622	12/2015	55.6 – 83.4	69.5

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27621	163.6	PASS
27622	65.0	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY: _____


Rokana Magdaleno
Quality Control

DATE: _____

07/24/2015