



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **TRIGLYCERIDE – GPO**      MANUFACTURE DATE: 11/18/2010  
 CATALOG #: T531-400  
 LOT #: 29098      EXP. DATE: 12/2012

COMPONENTS:    1. REAGENT      LOT#: 27058    EXP. DATE: 08/2013  
                          2. STANDARD      LOT#: 23965    EXP. DATE: 12/2012

BLANK:            0.018            AT 550 NM

STANDARD:      200 mg/dL

**TOLERANCE LEVELS OF SERUM CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> <u>(mg/dL)</u>
15271	03/2012	137 – 207	172
15272	03/2012	59 - 95	77

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	160.6	194.6	177.5	2/2
15272	86.5	78.5	82.2	2/2

COMMENT: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	<u>11/18/2010</u>
-----------------	--	-------	-------------------