



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URIC ACID** MANUFACTURE DATE: 11/3/2009
 CATALOG #: U580-240
 LOT #: 21307 EXP. DATE: 9/2011

<u>COMPONENTS</u>	<u>CATALOG #</u>	<u>LOT #</u>	<u>EXP. DATE</u>
1. URIC ACID REAGENT	U580-A	16815	9/2011
2. URIC ACID STANDARD (5.0 mg/dl)	U580-B	19147	4/2012

BLANK: 0.009 AT 505 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mg/dl)
15271	3/2012	4.2 – 6.2	5.2
15272	3/2012	6.9 – 10.3	8.7

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	4.3	5.7	5.0	2/2
15272	8.3	8.8	8.6	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Gascón* DATE: 11/3/2009
 Quality Control