

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URIC ACID** MANUFACTURE DATE: 08/18/2010
CATALOG #: U580-240
LOT #: 27127 EXP. DATE: 09/2011

<u>COMPONENTS</u>	<u>CATALOG #</u>	<u>LOT #</u>	<u>EXP. DATE</u>
1. URIC ACID REAGENT	U580-A	16815	09/2011
2. URIC ACID STANDARD (5.0 mg/dl)	U580-B	24717	03/2012

BLANK: 0.017 AT 505 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (mg/dl)</u>
15271	03/2012	3.9 - 6.1	5.0
15272	03/2012	6.6 - 10.0	8.3

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	5.4	5.1	5.3	2/2
15272	8.5	8.1	8.3	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Garcia* DATE: 08/18/2010
Quality Control