



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URIC ACID** MANUFACTURE DATE: 05/14/2014
 CATALOG #: U580-240
 LOT #: 54768 EXP. DATE: 01/2017

COMPONENTS: 1. URIC ACID REAGENT LOT #: 54188 EXP. DATE: 04/2017
 2. URIC ACID STANDARD LOT #: 52744 EXP. DATE: 01/2017

TOLERANCE LEVELS OF SERUM CONTROL


<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dl)
27611	05/2015	4.7 – 6.6	5.6
27602	03/2015	7.3 – 10.3	8.8

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27611	4.8	PASS
27602	8.7	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY:


 Jacqueline Sanchez
 Quality Control

DATE:

05/15/2014