



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **URS-10**  
 CATALOG#: **URS-10-100**  
 LOT #: **22556**

MANUFACTURE DATE: **11/24/2009**  
 EXP. DATE: **04/2012**

<u>COMPONENTS:</u>	<u>CATALOG#:</u>	<u>LOT #:</u>	<u>EXP. DATE:</u>
LEUKOCYTES	URSP-LEU	21452	05/2012
NITRITE	URSP-NIT	22461	05/2012
UROBILINOGEN	URSP-URO	21263	04/2012
PROTEIN	URSP-PRO	21530	05/2012
pH	URSP-PH	22131	05/2012
BLOOD	URSP-BLO	22360	05/2012
SPECIFIC GRAVITY	URSP-SG	21430	05/2012
KETONE	URSP-KET	22028	05/2012
BILIRUBIN	URSP-BIL	21147	04/2012
GLUCOSE	URSP-GLU	22130	05/2012

**TOLERANCE LEVELS OF CONTROLS**

Controls	<u>Leuko- cytes</u>	<u>Nitrite</u>	<u>Urobili- nogen</u>	<u>Protein</u>	<u>pH</u>	<u>Blood</u>	<u>Specific Gravity</u>	<u>Ketone</u>	<u>Bilirubin</u>	<u>Glucose</u>
KOVA I LOT#: 134042 EXP: 08/2010	2+ to 3+	POS (+)	4.0 to ≥8.0 EU/dl	100 to ≥300 mg/dl	7.0 to 8.0	2+ to 3+	1.010 to ≥1.030	40 to ≥80 mg/dl	2+ to 3+	500 to ≥1000 mg/dl
KOVA II LOT#: 133585 EXP: 05/2010	Trace to 2+	POS (+)	0.2 to 1.0 EU/dl	NEG (-) to 30 mg/dl	6.5 to 8.5	Trace to 2+	≤1.005 to 1.015	Trace to ≥40 mg/dl	1+ to 2+	Trace to 500 mg/dl
KOVA III LOT#: 134047 EXP: 05/2011	NEG (-)	NEG (-)	0.2 to 1.0 EU/dl	NEG (-)	5.0 to 6.5	NEG (-)	≤1.005 to 1.015	NEG (-)	NEG (-)	NEG (-)

**TEST VALUES/RESULTS**

<u>Controls</u>	<u>Leuko- cytes</u>	<u>Nitrite</u>	<u>Urobili- nogen</u>	<u>Protein</u>	<u>pH</u>	<u>Blood</u>	<u>Specific Gravity</u>	<u>Ketone</u>	<u>Bilirubin</u>	<u>Glucose</u>
KOVA I	2+ 2+	+ +	8.0 8.0	300 300	7.0 7.5	3+ 3+	1.015 1.015	80 80	3+ 3+	1000 1000
KOVA II	2+ 2+	+ +	0.2 0.2	15 15	7.5 7.5	1+ 1+	1.010 1.010	40 40	2+ 1+	250 250
KOVA III	- -	- -	0.2 0.2	- -	6.5 6.5	- -	1.005 1.005	- -	- -	- -

COMMENTS: VISUAL CHECK OK.

APPROVE  
D BY:

  
Quality Control

DATE 11/24/2009

: