



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **URS-10 100 Teco Finished Kits** MANUFACTURE DATE: 05/21/2010  
 CATALOG#: URS-10-100  
 LOT #: 25872 EXP. DATE: 08/2012

<u>COMPONENTS:</u>	<u>CATALOG#:</u>	<u>LOT #:</u>	<u>EXP. DATE:</u>
LEUKOCYTES	URSP-LEU	25537	09/2012
NITRITE	URSP-NIT	25538	09/2012
UROBILINOGEN	URSP-URO	25689	09/2012
PROTEIN	URSP-PRO	25714	09/2012
pH	URSP-PH	25743	09/2012
BLOOD	URSP-BLO	25645	09/2012
SPECIFIC GRAVITY	URSP-SG	25742	09/2012
KETONE	URSP-KET	25672	09/2012
BILIRUBIN	URSP-BIL	25608	09/2012
GLUCOSE	URSP-GLU	25253	08/2012

**TOLERANCE LEVELS OF CONTROLS**

Controls	<u>Leuko-cytes</u>	<u>Nitrite</u>	<u>Urobili-nogen</u>	<u>Protein</u>	<u>pH</u>	<u>Blood</u>	<u>Specific Gravity</u>	<u>Ketone</u>	<u>Bilirubin</u>	<u>Glucose</u>
KOVA I LOT#: 134042 EXP: 08/2010	2+ to 3+	POS (+)	4.0 to ≥8.0 EU/dl	100 to ≥300 mg/dl	7.0 to 8.0	2+ to 3+	1.010 to ≥1.030	40 to ≥80 mg/dl	2+ to 3+	500 to ≥1000 mg/dl
KOVA II LOT#: 136668 EXP: 02/2011	Trace to 2+	POS (+)	0.2 to 1.0 EU/dl	NEG (-) to 30 mg/dl	6.5 to 8.5	Trace to 2+	≤1.005 to 1.015	Trace to ≥40 mg/dl	1+ to 2+	Trace to 500 mg/dl
KOVA III LOT#: 137752 EXP: 02/2012	NEG (-)	NEG (-)	0.2 to 1.0 EU/dl	NEG (-)	5.0 to 6.5	NEG (-)	≤1.005 to 1.015	NEG (-)	NEG (-)	NEG (-)

**TEST VALUES/RESULTS**

<u>Controls</u>	<u>Leuko-cytes</u>	<u>Nitrite</u>	<u>Urobili-nogen</u>	<u>Protein</u>	<u>pH</u>	<u>Blood</u>	<u>Specific Gravity</u>	<u>Ketone</u>	<u>Bilirubin</u>	<u>Glucose</u>
KOVA I	3+ 3+	+ +	8.0 8.0	100 300	7.0 7.5	3+ 3+	1.010 1.015	80 80	3+ 3+	1000 1000
KOVA II	2+ 2+	+ +	0.2 0.2	- Trace	7.0 7.5	1+ 2+	1.010 1.015	40 40	2+ 2+	250 250
KOVA III	- -	- -	0.2 0.2	- -	6.0 5.0	- -	1.010 1.015	- -	- -	- -

COMMENTS: VISUAL CHECK OK.

APPROVED  
 BY:

*V. Arseny*  
 Quality Control

DATE: 05/21/2010