



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URS-10 100 Teco Finished Kit** MANUFACTURED DATE: 08/24/2011
 CATALOG#: URS-10-100
 LOT #: 34561 EXP. DATE: 08/2013

COMPONENTS:	CATALOG#:	LOT #:	EXP. DATE:
LEUKOCYTES	URSP-LEU	33072	10/2013
NITRITE	URSP-NIT	32006	08/2013
UROBILINOGEN	URSP-URO	32191	08/2013
PROTEIN	URSP-PRO	31765	08/2013
pH	URSP-PH	31829	08/2013
BLOOD	URSP-BLO	32169	08/2013
SPECIFIC GRAVITY	URSP-SG	31848	08/2013
KETONE	URSP-KET	31824	08/2013
BILIRUBIN	URSP-BIL	32007	08/2013
GLUCOSE	URSP-GLU	31974	08/2013

TOLERANCE LEVELS OF CONTROLS

Controls	Leuko-cytes	Nitrite	Urobili-nogen	Protein	pH	Blood	Specific Gravity	Ketone	Bilirubin	Glucose
KOVA I LOT#:140399-8 EXP:06/2012	2+ to 3+	POS (+)	4.0 to ≥8.0 EU/dl	100 to ≥2000 mg/dl	7.0 to 8.5	2+ to 3+	1.010 to ≥1.030	40 to ≥160 mg/dl	2+ to 3+	500 to ≥2000 mg/dl
KOVA II LOT#:139317-7 EXP: 02/2012	Trace to 2+	POS (+)	0.2 to 1.0 EU/dl	NEG (-) to 30 mg/dl	6.5 to 8.5	Trace to 2+	≤1.005 to 1.015	Trace to ≥40 mg/dl	1+ to 2+	Trace to 500 mg/dl
KOVA III LOT#:141287-1 EXP: 08/2013	NEG (-)	NEG (-)	0.2 to 1.0 EU/dl	NEG (-)	5.0 to 6.5	NEG (-)	≤1.005 to 1.015	NEG (-)	NEG (-)	NEG (-)

TEST VALUES/RESULTS

Controls	Leuko-cytes	Nitrite	Urobili-nogen	Protein	pH	Blood	Specific Gravity	Ketone	Bilirubin	Glucose
KOVA I	2+ 2+	+ +	8.0 8.0	300 300	8.0 7.5	3+ 3+	1.015 1.015	40 80	3+ 3+	1000 1000
KOVA II	1+ 2+	+ +	0.2 0.2	- -	8.0 7.5	Trace 2+	1.010 1.010	15 Trace	1+ 2+	100 100
KOVA III	- -	- -	0.2 0.2	- -	6.5 5.0	- -	1.010 1.010	- -	- -	- -

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:08/25/2011
--------------	--	-----------------