



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URS-10 CLICHEK F/KIT**
 CATALOG#: **URS-10-100CK**
 LOT #: **27786**

MANUFACTURE DATE: **09/30/2010**

EXP. DATE: **07/2012**

COMPONENTS:	LEUKOCYTES	LOT #:	26151	EXP. DATE:	10/2012
	NITRITE	LOT #:	26153	EXP. DATE:	10/2012
	UROBILINOGEN	LOT #:	26181	EXP. DATE:	10/2012
	PROTEIN	LOT #:	26049	EXP. DATE:	10/2012
	pH	LOT #:	26237	EXP. DATE:	10/2012
	BLOOD	LOT #:	26003	EXP. DATE:	10/2012
	SPECIFIC GRAVITY	LOT #:	26221	EXP. DATE:	10/2012
	KETONE	LOT #:	24902	EXP. DATE:	07/2012
	BILIRUBIN	LOT #:	26159	EXP. DATE:	10/2012
	GLUCOSE	LOT #:	24874	EXP. DATE:	07/2012

TOLERANCE LEVELS OF CONTROLS

Controls	Leuko-cytes	Nitrite	Urobili-nogen	Protein	pH	Blood	Specific Gravity	Ketone	Bilirubin
KOVA I LOT#: 139036-4 EXP:01/2012	2+ to 3+	POS (+)	4.0 to ≥8.0 EU/dl	100 to ≥300 mg/dl	7.0 to 8.0	2+ to 3+	1.010 to ≥1.030	40 to ≥80 mg/dl	2+ to 3+
KOVA II LOT#:138183-2 EXP: 07/2011	Trace to 2+	POS (+)	0.2 to 1.0 EU/dl	NEG (-) to 30 mg/dl	6.5 to 8.5	Trace to 2+	≤1.005 to 1.015	Trace to ≥40 mg/dl	1+ to 2+
KOVA III LOT#:137752-3 EXP: 02/2012	NEG (-)	NEG (-)	0.2 to 1.0 EU/dl	NEG (-)	5.0 to 6.5	NEG (-)	≤1.005 to 1.015	NEG (-)	NEG (-)

Glucose:

Control Level:	Level 1	Level 2	Level 3	Level4	Level 5
Target (mg/dL):	0	50	100	300	≥1000
Range (mg/dL):	0	25 - 75	50 - 150	150 - 400	500 - ≥1000


TEST VALUES/RESULTS

Controls	Leukocytes	Nitrite	Urobili-nogen	Protein	pH	Blood	Specific Gravity	Ketone	Bilirubin
KOVA I	2+ 2+	+ +	8.0 8.0	300 300	7.5 7.5	2+ 2+	1.015 1.015	80 80	3+ 3+
KOVA II	1+ 1+	+ +	0.2 0.2	Trace Trace	8.0 8.0	1+ 1+	1.010 1.010	40 40	2+ 2+
KOVA III	- -	- -	0.2 0.2	- -	6.5 6.5	- -	1.005 1.005	- -	- -

Glucose:

Test Results:	0	50	100	300	1000
---------------	---	----	-----	-----	------

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	<u>09/30/2010</u>
--------------	--	-------	-------------------