



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **URS-10 100 CLICHEK Fin. Kit**

MANUFACTURE DATE: 03/15/2011

CATALOG#: URS-10-100CK

LOT #: 31295

EXP. DATE: 05/2013

COMPONENTS:	LEUKOCYTES	LOT #:	30104	EXP. DATE:	05/2013
	NITRITE	LOT #:	30448	EXP. DATE:	06/2013
	UROBILINOGEN	LOT #:	29876	EXP. DATE:	05/2013
	PROTEIN	LOT #:	30450	EXP. DATE:	06/2013
	pH	LOT #:	29899	EXP. DATE:	05/2013
	BLOOD	LOT #:	30229	EXP. DATE:	05/2013
	SPECIFIC GRAVITY	LOT #:	30010	EXP. DATE:	05/2013
	R.KETONE	LOT #:	30503	EXP. DATE:	06/2013
	BILIRUBIN	LOT #:	30584	EXP. DATE:	06/2013
	R.GLUCOSE	LOT #:	29290	EXP. DATE:	05/2013

**TOLERANCE LEVELS OF CONTROLS**

Controls	Leuko-cytes	Nitrite	Urobili-nogen	Protein	pH	Blood	Specific Gravity	Ketone	Bilirubin
KOVA I LOT#: 139036-4 EXP:01/2012	2+ to 3+	POS (+)	4.0 to ≥8.0 EU/dl	100 to ≥2000 mg/dl	7.0 to 8.5	2+ to 3+	1.010 to ≥1.030	40 to ≥160 mg/dl	2+ to 3+
KOVA II LOT#:139317-7 EXP: 02/2012	Trace to 2+	POS (+)	0.2 to 1.0 EU/dl	NEG (-) to 30 mg/dl	6.5 to 8.5	Trace to 2+	≤1.005 to 1.015	Trace to ≥40 mg/dl	1+ to 2+
KOVA III LOT#:140400-1 EXP: 03/2013	NEG (-)	NEG (-)	0.2 to 1.0 EU/dl	NEG (-)	5.0 to 6.5	NEG (-)	≤1.005 to 1.015	NEG (-)	NEG (-)

Glucose:

Control Level:	Level 1	Level 2	Level 3	Level4	Level 5
Target (mg/dL):	0	50	100	300	≥1000
Range (mg/dL):	0	25 - 75	50 - 150	150 - 400	500 ->1000

**TEST VALUES/RESULTS**

Controls	Leukocytes	Nitrite	Urobili-nogen	Protein	pH	Blood	Specific Gravity	Ketone	Bilirubin
KOVA I	3+ 3+	+ +	8.0 8.0	300 300	8.0 8.0	3+ 3+	1.020 1.020	80 80	3+ 3+
KOVA II	1+ 1+	+ +	0.2 0.2	Trace Trace	7.5 7.5	1+ 1+	1.010 1.010	40 40	1+ 1+
KOVA III	- -	- -	0.2 0.2	- -	6.5 6.5	- -	1.005 1.005	- -	- -

Glucose:

Test Results:	0	50	100	300	1000
---------------	---	----	-----	-----	------

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	03/15/2011
--------------	--	-------	------------