



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URS-10 100 CLICHEK**

MANUFACTURE DATE: 07/29/2011

Finished Kit

CATALOG#: URS-10-100CK

LOT #: 33866

EXP. DATE: 07/2013

COMPONENTS:	LEUKOCYTES	LOT #:	33625	EXP. DATE:	11/2013
	NITRITE	LOT #:	33748	EXP. DATE:	11/2013
	UROBILINOGEN	LOT #:	32320	EXP. DATE:	09/2013
	PROTEIN	LOT #:	33294	EXP. DATE:	10/2013
	pH	LOT #:	32316	EXP. DATE:	09/2013
	BLOOD	LOT #:	32169	EXP. DATE:	08/2013
	SPECIFIC GRAVITY	LOT #:	32319	EXP. DATE:	09/2013
	R. KETONE	LOT #:	31324	EXP. DATE:	07/2013
	BILIRUBIN	LOT #:	32007	EXP. DATE:	08/2013
	R. GLUCOSE	LOT #:	31325	EXP. DATE:	07/2013

TOLERANCE LEVELS OF CONTROLS

Controls	<u>Leuko-</u> <u>cytes</u>	<u>Nitrite</u>	<u>Urobili-</u> <u>nogen</u>	<u>Protein</u>	<u>pH</u>	<u>Blood</u>	<u>Specific</u> <u>Gravity</u>	<u>Ketone</u>	<u>Bilirubin</u>
KOVA I LOT#: 134042 EXP:08/2010	2+ to 3+	POS (+)	4.0 to ≥8.0 EU/dl	100 to ≥2000 mg/dl	7.0 to 8.5	2+ to 3+	1.010 to ≥1.030	40 to ≥160 mg/dl	2+ to 3+
KOVA II LOT#:133585 EXP: 05/2010	Trace to 2+	POS (+)	0.2 to 1.0 EU/dl	NEG (-) to 30 mg/dl	6.5 to 8.5	Trace to 2+	≤1.005 to 1.015	Trace to ≥40 mg/dl	1+ to 2+
KOVA III LOT#:134047 EXP: 05/2011	NEG (-)	NEG (-)	0.2 to 1.0 EU/dl	NEG (-)	5.0 to 6.5	NEG (-)	≤1.005 to 1.015	NEG (-)	NEG (-)

Glucose:

Control Level:	Level 1	Level 2	Level 3	Level4	Level 5
Target (mg/dL):	0	50	100	300	≥1000
Range (mg/dL):	0	25 - 75	50 - 150	150 - 400	500 - ≥1000

TEST VALUES/RESULTS

<u>Controls</u>	<u>Leukocytes</u>	<u>Nitrite</u>	<u>Urobili-</u> <u>nogen</u>	<u>Protein</u>	<u>pH</u>	<u>Blood</u>	<u>Specific</u> <u>Gravity</u>	<u>Ketone</u>	<u>Bilirubin</u>
KOVA I	2+ 2+	+ +	8.0 8.0	300 300	7.5 7.5	3+ 3+	1.020 1.020	80 80	3+ 3+
KOVA II	1+ 1+	+ +	0.2 0.2	Trace Trace	7.0 7.0	2+ 2+	1.010 1.010	40 40	2+ 2+
KOVA III	- -	- -	0.2 0.2	- -	6.0 5.0	- -	1.010 1.010	- -	- -

Glucose:

Test Results:	0	50	100	300	≥1000
---------------	---	----	-----	-----	-------

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	07/29/2011
-----------------	--	-------	------------