



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URS-1G 100 strips** MANUFACTURE DATE: 1/4/10
Teco Finished Kit
 CATALOG#: **URS-1G-100**
 LOT #: 23288 EXP. DATE: 5/12

COMPONENT: URS PAPER - GLUCOSE LOT #: 22130 EXP. DATE: 05/2012


TOLERANCE LEVELS OF CONTROLS

CONTROLS	LOT #	EXP. DATE	GLUCOSE RANGE
KOVA I	134042	08/2010	500 to \geq 1000 mg/dl
KOVA II	133585	05/2010	Trace to 500 mg/dl
KOVA III	134047	05/2011	Negative (-)

TEST VALUES/RESULTS

CONTROLS	RESULT 1	RESULT 2	# PASS
LEVEL I	1000	1000	2 / 2
LEVEL II	250	250	2 / 2
LEVEL III	-	-	2 / 2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	<u>1/4/10</u>
-----------------	--	-------	---------------