

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URS-1G 100 strips** MANUFACTURE DATE: 08/06/2010
 Teco Finished Kit
CATALOG#: **URS-1G-100**
LOT #: 26983 EXP. DATE: 10/2012

COMPONENT: URS PAPER - GLUCOSE LOT #: 26152 EXP. DATE: 10/2012


TOLERANCE LEVELS OF CONTROLS

CONTROLS	LOT #	EXP. DATE	GLUCOSE RANGE
KOVA I	134042	08/2010	500 to ≥1000 mg/dl
KOVA II	136668	02/2011	Trace to 500 mg/dl
KOVA III	137752	02/2012	Negative (-)

TEST VALUES/RESULTS

CONTROLS	RESULT 1	RESULT 2	# PASS
LEVEL I	1000	1000	2 / 2
LEVEL II	250	250	2 / 2
LEVEL III	-	-	2 / 2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:08/06/2010	
--------------	--	-----------------	--