

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(BATCH RELEASE)

PRODUCT NAME: **URS-1P 100 bottle / No label** MANUFACTURE DATE: 09/16/2010
CATALOG #: URS-1P-100B
LOT #: 27658 EXP. DATE: 01/2013

COMPONENT: URS Paper - Protein LOT #: 27552 EXP. DATE: 01/2013

TOLERANCE LEVELS OF CONTROLS

CONTROLS	LOT #	EXP. DATE	PROTEIN
Kova-Trol I	139036-4	01/2012	100 to \geq 300 mg/dL
Kova-Trol II	138183-2	07/2011	Negative (-) to 30 mg/dL
Kova-Trol III	137752-3	02/2012	Negative (-)

TEST VALUES/RESULTS

CONTROLS	1	2	# PASS
LEVEL I	300	300	2 / 2
LEVEL II	Trace	Trace	2 / 2
LEVEL III	-	-	2 / 2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 09/16/2010
Quality Control