



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **URS-2P 100 Teco Finished Kit** MANUFACTURE DATE: 10/26/2010  
 CATALOG #: **URS-2P-100**  
 LOT #: 28586 EXP. DATE: 02/2013

COMPONENTS: URS PAPER – GLUCOSE LOT #: 28118 EXP. DATE: 02/2013  
 URS PAPER – PROTEIN LOT #: 28117 EXP. DATE: 02/2013

**TOLERANCE LEVELS OF CONTROLS**

CONTROLS	LOT #	EXP. DATE	GLUCOSE	PROTEIN
KOVA I	139036-4	01/2012	500 to $\geq$ 1000 mg/dl	100 to $\geq$ 300 mg/dl
KOVA II	138183-2	07/2011	Trace to 500 mg/dl	Negative (-) to 30 mg/dl
KOVA III	139316-9	08/2012	Negative (-)	Negative (-)

**TEST VALUES/RESULTS**

CONTROLS	GLUCOSE		PROTEIN		# PASS
	1	2	1	2	
KOVA I	500	500	300	300	4/4
KOVA II	250	250	Trace	Trace	4/4
KOVA III	-	-	-	-	4/4

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 10/27/2010  
 Quality Control