



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URS-3**
CATALOG #: URS-3-100
LOT #: 24414

MANUFACTURE DATE: 03/10/2010
EXP. DATE: 07/2012

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT #:</u>	<u>EXP. DATE:</u>
URS Paper - Protein	URSP-PRO	24031	08/2012
URS Paper - Glucose	URSP-GLU	24030	08/2012
URS Paper - pH	URSP-PH	22781	07/2012

TOLERANCE LEVELS OF CONTROLS


CONTROLS	LOT #	EXP. DATE	PROTEIN	GLUCOSE	pH
KOVA I	134042	08/2010	100 to $\geq$ 300 mg/dl	500 to $\geq$ 1000 mg/dl	7.0 – 8.0
KOVA II	133585	05/2010	Negative (-) to 30 mg/dl	Trace to 500 mg/dl	6.5 – 8.5
KOVA III	134047	05/2011	Negative (-)	Negative (-)	5.0 – 6.5

TEST VALUES/RESULTS

CONTROLS	PROTEIN		GLUCOSE		pH		# PASS
	<u>Result 1</u>	<u>Result 2</u>	<u>Result 1</u>	<u>Result 2</u>	<u>Result 1</u>	<u>Result 2</u>	
LEVEL I	300	300	1000	1000	7.5	7.5	6 / 6
LEVEL II	Trace	Trace	250	250	8.0	8.0	6 / 6
LEVEL III	-	-	-	-	6.5	6.5	6 / 6

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:


Quality Control

DATE: 03/10/2010